
Instructions for Use | Gebrauchsanleitung | Manuel d'utilisation | Istruzioni per l'uso | Instrucciones de uso | Gebruiksaanwijzing

Clipster®

Aseptic Disconnector | Sterile Trennvorrichtung | Déconnecteur aseptique | Disconnettore asettico | Desconectador aséptico | Steriele scheidingsinrichting

85032-538-41

SARTORIUS

English	page	3
Deutsch	Seite	29
Français	page	56
Italiano	pagina	84
Español	página	111
Nederlands	Pagina	138

Contents

1 About These Instructions 5
1.1 Validity 5
1.2 Target Groups 5
1.3 Symbols Used 6
1.3.1 Warnings in Operation Descriptions 6
1.3.2 Other Symbols Used 6

2 Safety Instructions 6
2.1 Intended Use 6
2.1.1 Modifications to the Product 7
2.1.2 Repairs and Maintenance on the Product 7
2.2 Qualifications of Personnel 8
2.3 Functionality of the Product Parts 8
2.4 Safety Equipment 8
2.5 Safety Information on the Product 8
2.6 Consumables and Spare Parts 9
2.7 Personal Protective Equipment 9
2.8 Heavy Product Weight 9
2.9 Escaping Liquids 9

3 Product Description 10
3.1 Product Overview 10
3.2 Overview of the Clipster® Aseptic Disconnector 11
3.3 Safety Equipment 12
3.4 Product Symbols 12

4 Operation 13
4.1 Tips for Pre-assembling Tubing 13
4.2 Disconnecting Tubing 14

5 Cleaning and Maintenance 20
5.1 Cleaning 20
5.2 Maintenance Schedule 20
5.3 Replacing the Blade 21

6 Disposal 23
6.1 Decontaminating 23
6.2 Disposing of Product and Parts 23

7 Technical Data **24**
7.1 Dimensions and Weights 24
7.2 Materials 24
7.3 Cleaning Agents and Cleaning Procedures 24
7.4 Storage Conditions 25

8 Consumables and Spare Parts **26**
8.1 Consumables 26
8.2 Spare Parts 26

9 Sartorius Service **28**

1 About These Instructions

1.1 Validity

These instructions are part of the product, they must be read in full and stored. These instructions apply to the following versions of the product:

Hand Held Tool	Material	Item Number
Clipster® XS \| S		642-S--T
Clipster® L		642-L--T

Clipster® Aseptic Disconnector	Material	Item Number
Size XS	Polycarbonate	642XS-------D
		642-XS--200---N
Size S	Polycarbonate	642-S-------D
		642-S--200---N
Size L	Polycarbonate	642-L-------D
		642-L--200---N

1.2 Target Groups

These instructions are addressed to the following target groups. The target groups must possess the knowledge specified below.

Target Group	Knowledge and Qualifications
User	The user is familiar with the product and the associated work processes. The user is familiar with the hazards that may arise when working with the product, and is able to avoid these hazards.

1.3 Symbols Used

1.3.1 Warnings in Operation Descriptions

⚠ CAUTION
Denotes a hazard that may result in moderate or minor injury if it is not avoided.

1.3.2 Other Symbols Used

▶ Required action: Describes activities that must be carried out. The activities in the sequence must be carried out in succession.

▷ Result: Describes the result of the activities carried out.

Figures in These Instructions

Depending on the product configuration, the figures may differ slightly from the supplied product. The variants shown in these instructions are examples.

2 Safety Instructions

2.1 Intended Use

The Clipster® Disconnector is intended for sterile disconnection of tubing in biopharmaceutical applications.

The product is intended exclusively for use in accordance with these instructions. Any other use is considered improper.

Foreseeable Misuse

The following applications are **n**ot permitted:

- Disconnection of unauthorized tubing
- Disconnection of unsuitable tubing, e.g., braided tubing
- Disconnection of tubing without using the Clipster® Disconnector

Operating Conditions for the Product

Only use the product indoors and under usual laboratory conditions.

2.1.1 Modifications to the Product

If the product is modified: Persons may be put at risk. Product-specific documents and product approvals may lose their validity.

For queries regarding modifications to the product, contact Sartorius.

2.1.2 Repairs and Maintenance on the Product

Product repairs and maintenance may only be carried out by persons with specialized knowledge of the product. If the product is not repaired or maintained by a specialist: Persons may be put at risk. Product-specific documents and product approvals may lose their validity.

Sartorius recommends that any repair work, even that carried out after the end of the warranty period, is carried out by Sartorius Service or after consulting with Sartorius Service.

Only the maintenance tasks described in these instructions should be carried out. For maintenance tasks that need to be carried out by Sartorius Service, contact Sartorius Service.

2.2 Qualifications of Personnel

Persons who do not possess adequate knowledge about how to use the product safely may injure themselves and other persons.

If a particular qualification is required for an activity: The target group will be specified. If no qualification is specified: The activity may be carried out by the “user” target group.

2.3 Functionality of the Product Parts

Non-functioning product parts, e.g., as a result of damage or wear, can cause malfunctions. There is a risk of injury to persons.

- ▶ If product parts are not functioning: Do not use the product.
- ▶ Comply with the maintenance intervals (intervals and maintenance tasks see Chapter “5.2 Maintenance Schedule”, page 20).

2.4 Safety Equipment

The product’s safety equipment protects persons who work with the product against the hazards associated with it. If the product’s safety equipment is detached or modified: People may be seriously injured.

- ▶ Do not dismantle, modify, or disable the safety equipment (for safety equipment, see Chapter “3.3 Safety Equipment”, page 12).

2.5 Safety Information on the Product

Symbols, e.g., warnings and safety stickers, are safety information for handling the product. Missing or illegible safety information may result in this information not being observed. There is a risk of injury to persons.

- ▶ Do not conceal, remove, or modify the symbols.
- ▶ Replace the symbols if they become illegible.

2.6 Consumables and Spare Parts

The use of unsuitable consumables and spare parts can affect the functionality and safety of the product and have the following consequences:

- Risk of injury to persons
- Damage, malfunctions, or failure of the product

▶ Only use consumables and spare parts that have been approved by Sartorius and are safe for use.

2.7 Personal Protective Equipment

Personal protective equipment protects against hazards arising from the product. If the personal protective equipment is missing or is unsuitable for the work processes on the product: Persons may be injured.

The following personal protective equipment must be worn:

- Protective work clothing
- Safety gloves
- Safety glasses
- Safety boots

2.8 Heavy Product Weight

The product is very heavy. When lifting, transporting, and operating the product, there is a risk of injury, e.g., as a result of the product falling.

▶ Wear the necessary personal protective equipment.

2.9 Escaping Liquids

If tubing other than the approved tubing or tubing materials is used: This may result in the tubing not separating properly. Liquids may leak from the tubing.

▶ Only use approved tubing and tubing materials.
▶ Depending on the liquids used: Wear the necessary personal protective equipment.

3 Product Description

3.1 Product Overview


Fig. 1: Overview of hand held tool for Clipster® L

Pos.	Name	Description
1	Cutting handle	Used to operate the blade.
2	Blade	Separates the Clipster® Disconnector.
3	Clipster® seat	Holds the Clipster® Disconnector.
4	Pressure plate	Presses the Clipster® Disconnector firmly into the Clipster® seat.
5	Clamping handle	Used to operate the clamp.
6	Body	
7	Tool handle	Used to hold the product.

3.2 Overview of the Clipster® Aseptic Disconnector


Fig. 2: Components of the Clipster® Aseptic Disconnector

Pos.	Description
1	Bottom part of the Clipster® Aseptic Disconnector
2	Top part of the Clipster® Aseptic Disconnector


Fig. 3: Latching positions of the Clipster® Disconnector

Pos.	Description
1	Latching position 1
2	Latching position 2

3.3 Safety Equipment


Fig. 4: Safety equipment

Pos.	Name	Description
1	Protective cover	Protects against direct contact with the blade.
2	Blade safety switch	Used to lock the blade.

3.4 Product Symbols

Symbol	Meaning
	Hand injuries: Risk of cuts from the sharp blade


4 Operation

4.1 Tips for Pre-assembling Tubing

Tubing is not always perfectly round and may be oval in shape. If oval tubing is used, the following procedure is recommended.

Procedure


▶ Place the oval tubing into the bottom part of the Clipster® Aseptic Disconnector in such a way that the tubing lies flat.


▷ This makes it easier to then snap the top part of the Clipster® Aseptic Disconnector to the bottom part in the center.

4.2 Disconnecting Tubing

Requirements


– The cutting handle (1) and the clamping handle (2) should be pointing upwards.


– The blade safety switch (1) must be engaged as shown.


Procedure

▶ Slide together the two bottom parts of the Clipster® Disconnector along the grooved surfaces.



▶ Place the tubing in the two bottom parts.

▶ Snap the first top part into latching position 1 of the bottom part.

▶ Check whether the top part is centered. If the top part is not centered: Center the tubing as shown.

▶ Snap the second top part into latching position 1 of the second bottom part.
▶ Visually check that the Sartorius logo on one top part is aligned with the arrow on the other top part.
▶ Check whether the second top part is centered. If the top part is not centered: Center the tubing.

Only necessary for Tuflux SIL and TPE 1/2″ x 3/4″ tubing


To make it easier to decouple Tuflux SIL and TPE 1/2" x 3/4" tubing, the friction between the tubing and the Clipster® Disconnector can be reduced.

Please ensure to use the new version of the Clipster® L-Tool (lot number starting from 4222).


Procedure

- ▶ Carefully stretch the tubing before inserting the Clipster® Disconnector into the hand held tool.
 - ▶ To do this, hold the Clipster® Disconnector as shown. Place your thumbs under the bottom parts and your index fingers on the side of the top parts.

- ▶ Stretch the tubing by performing a bending motion with your hands.
- ▶ Keeping your thumbs in the same position below the bottom parts, move your index fingers away from the top parts until...

- ▶ ...there is a minimum distance of 1 cm between the top parts and your index fingers.

Disconnecting Tubing


⚠ CAUTION

Risk of leakage when disconnecting pressurized tubing!


Disconnection can be performed on empty and partly filled tubing.

▶ When disconnecting partly filled tubing, ensure no pressure is generated on the tubing.


Procedure


▶ Insert the assembled Clipster® Disconnector top side (with the Sartorius logo) first into the Clipster® seat.


▷ The Clipster® Disconnector must be positioned in the Clipster® seat as shown.


- ▶ Pull the clamping handle downwards to clamp the Clipster® Disconnector securely.
- ▷ An audible click indicates that the Clipster® Disconnector is seated correctly.


- ▶ Turn the blade safety switch to the left to release the blade.

- ▶ Hold the clamping handle facing downwards.
- ▶ While doing so, pull the cutting handle downwards to cut the tubing.


- ▶ Return the cutting handle to the starting position by pulling it upwards.
- ▶ Turn the blade safety switch back to the secured position.


▶ Pull the clamping handle back upwards to its starting position.

▶ Remove the Clipster® Disconnector from the Hand Held Tool.

▶ Push the two ends together (1) and slide them in opposite directions vertically against one another (2).

Only for sizes XS and S

▶ The Clipster® Disconnector can be separated by bending the two Clipster® parts.

▶ Check the two ends after disconnection.

⚠ CAUTION

Danger of Leakage

Where there is mechanical stress, e.g., pulling on the tubing, this may result in leaks and in liquid escaping.

▶ Do not subject the ends of the tubing to mechanical stress after disconnection.

5 Cleaning and Maintenance

5.1 Cleaning

The Clipster® Disconnector can be cleaned with most disinfection solutions such as isopropyl alcohol.

Procedure

▶ Only use suitable cleaning agents and cleaning procedures and observe the product information for the cleaning agent used (for cleaning agent, see Chapter 7.3, page 24).

▶ Wipe the housing with a slightly damp cloth. In the case of more heavy contamination, use a mild soapy solution or a suitable cleaning agent.

5.2 Maintenance Schedule

Interval	Component	Activity	Chapter, Page	Target Group
After approx. 100 disconnection operations	Blade	Replacing the blade	5.3, 21	User

5.3 Replacing the Blade


⚠ CAUTION

Risk of injury from the sharp blade!


When replacing the blade, there is direct contact with the blade.

- ▶ Handle the blade with extreme caution.
- ▶ Wear the necessary personal protective equipment.

Procedure


- ▶ Loosen the knurled screw (1) slightly.
- ▶ Lift the protective cover upwards (2).


- ▶ Unscrew the knurled screw (1) until the blade is released.
- ▶ Carefully pull out the blade (2).

- ▶ Insert a new blade.

▶ Tighten the knurled screw (1) to secure the blade.

▶ Close the protective cover (1).
▶ Tighten the knurled screw (2) on the protective cover.

6 Disposal

6.1 Decontaminating

Procedure

▶ Perform decontamination.

6.2 Disposing of Product and Parts

The product must be disposed of properly by disposal facilities.

Requirements

The product is decontaminated.

Procedure

▶ Dispose of the product in accordance with local government regulations.
▶ Dispose of the packaging | consumables in accordance with local government regulations.

7 Technical Data

7.1 Dimensions and Weights

	Unit	Value
Dimensions (W x L x H)	mm	58 x 262 x 235
Weight, approx.	kg	1.1

7.2 Materials

Material	Component
Stainless steel and plastic	Cutting handle, clamping handle
Stainless steel	Blade, pressure plate
Aluminium	Clipster® seat, tool handle, body
Polycarbonate	Protective cover

7.3 Cleaning Agents and Cleaning Procedures

Cleaning agents
– 70% ethanol \| 30% WFI – 70% isopropanol \| 30% WFI
No corrosive or scouring components
No chlorine or chloride-containing components
No solvents
Cleaning procedures
Wipe surfaces with soaked cloths

Remove residues with plastic cleaning brushes
No autoclaving
No spraying processes for cleaning or disinfection
No high-pressure cleaners or steam jet cleaners

7.4 Storage Conditions

Storage conditions of blades and spare parts
Store in original packaging, dry and at ambient temperature.
The blades and spare parts are made of stainless steel and plastic. Since blades and spare parts are not exposed to light and are never in contact with the liquid which remains inside the tubing, they have no expiry date.

8 Consumables and Spare Parts

8.1 Consumables

This table contains a selection of the consumables that can be ordered. For information on other products, please contact Sartorius.

Item	Order Number
Blades for Clipster® XS \| S tool	642-S--B
Blades for Clipster® L tool	642-L--B

8.2 Spare Parts


This table contains a selection of the spare parts that can be ordered. For information on other products, please contact Sartorius.

Item	Order Number
Hand held tool for Clipster® XS \| S	642-S--T
Hand held tool for Clipster® L	642-L--T
Spare parts box for hand held tool for Clipster® XS \| S	642-S-SP
Spare parts box for hand held tool for Clipster® L	642-L-SP

* Note: Tubing with a minimum length of 300 mm is required for the Clipster® Disconnector.

** C-Flex® is a registered trademark of Saint-Gobain Performance Plastics Corporation.

Scope of delivery of the spare parts box for hand held tool for Clipster®

Pos.	Item	Quantity
1	Plate for Clipster® seat	2
2	Cover	1
3	Screw for plate	4
4	Knurled screw for cover	1
5	Knurled screw for blade	1
6	Screw for cover	2
7	Sleeve for cover	2

9 Sartorius Service

Sartorius Service is available for queries regarding the product. Please visit the Sartorius website (www.sartorius.com) for information about the service addresses, services provided, or to contact a local representative.

For inquiries about the system or when contacting Sartorius Service in the event of a malfunction, ensure that you have the product information, e.g., serial number, close at hand so you can pass it on to Sartorius Service. Consult the information on the manufacturer's ID label.

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Phone: +49 551 308 0
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.
Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.
Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote all genders.

Copyright notice:
These instructions, including all components, are protected by copyright.
Any use beyond the limits of the copyright law is not permitted without our approval.
This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:

04 | 2026

© 2026
Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany

AM | Publication No.: SLO6006af260412

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** C-Flex® is een geregistreerd handelsmerk van Saint-Gobain Performance Plastics Corporation.